


pH-day 2



Double channel PHmeter
for gastroesophageal studies

Double channel PHmeter
for gastroesophageal studies

- User friendly.
- Small and light.
- Powerful

Menfis Division

Divisional Office:

Via della Beverara, 46/D - 40131 Bologna Italy.

Tel +39 051 635 10 74 - Fax +39 051 635 10 81

E-mail: info@menfis.it - Web site: www.menfis.it

MEDICA S.p.A.

HARDWARE DATA

- Duration of the recording: more than 24 hours;
 - High sampling rate at 1 sample/second;
 - Long duration memory even without batteries;
 - Postural sensor for automatic recording of the patient postural changes;
 - Three event marker - buttons;
 - Two LEDs for setting and reading;
 - Compatible with glass and antimony probes, single use and reusable;
 - Data transfer to the Personal Computer safe and quick by fiber-optical cable;
 - Range: 0 - 14 pH;
-
- Dimensions 10 x 6,5 x 2,5 cm;
 - Weight with batteries 140 g;
 - Powered by two alkaline batteries AA of 1,5 V;

SOFTWARE DATA

- Developed in Windows environment;
- Patient and studies database;
- Medical reports in Word, and in MS-OFFICE software packages;
- Several possibilities of analysis: Acid Reflux, Alkaline Reflux and Mixed Reflux with two channels;
- Calculation of the pHmeter parameters including the calculation of the DeMeester's score;
- Paediatric analysis with Boix-Ochoa reference table;
- Calculation of area under H+;
- Calculation of area under pH;
- Calculation of the Symptom Index and Symptom Sensitivity Index;
- Load curve of pH and H+;
- Oscillatoy Index;

** Technical data can be changed without any notice*

**Menfis
Division**

Divisional Office:
Via della Beverara, 46/D - 40131 Bologna Italy.
Tel +39 051 635 10 74 - Fax +39 051 635 10 81
E-mail: info@menfis.it - Web site: www.menfis.it

CE
0476

MEDICA S.p.A.

via degli artigiani, 7
41036 MEDOLLA MODENA - ITALY