


# VitaSmart

Simplicity. Value.

## The First Multi-Organ Pump for Hypothermic Oxygenated Machine Perfusion

- Your everyday liver and kidney pump
- Simple to operate
- Priced for routine use
- For pre-implant, in-hospital use
- CE marked
- Available in Europe from 2019



**Two uses** - able to pump liver or kidney

**Two essentials** - hypothermic and oxygenated to minimise damage

**Two oxygenators** - to allow either HOPE or DHOPE for livers

**Two controls** - select pressure and flow rates

## Disposable Kits

- Kidney
- Liver: HOPE Protocol (portal vein)
- Liver: DHOPE Protocol (portal vein and hepatic artery)





## Solutions

- Storage at room temperature (2-25°C)
- CE marked (Class III medical device)
- Stocks held in London, Paris and the Netherlands for rapid delivery



### Belzer UW<sup>®</sup> Cold Storage Solution

For flushing and preservation of abdominal organs (liver, kidney, pancreas)

- Eliminates Refrigeration Concerns
- No need to filter before use




### Belzer MPS<sup>®</sup> UW Machine Perfusion Solution

For machine perfusion and preservation of explanted kidneys

- Less viscous & lower potassium than Belzer UW<sup>®</sup> Cold Storage Solution
- For use in a variety of pumping machines
- Widely used in liver & kidney multi-national studies (COPE, HOPE, DHOPE)

## French-Sized Cannulae

- Suture bump for secure suturing
- Adjustable bump (larger cannulae only) to allow degree of flexibility of fit
- CE marked, colour-coded



Bridge to Life (Europe) Ltd. | For More Information Contact: Mike Kappler

Tel: +44(0)20 3411 8326 | Mobile: +44(0)77 6652 7191 | [www.bridgetolife.eu](http://www.bridgetolife.eu)

LU 311 The Light Bulb | 1 Filament Walk | London SW18 4GQ

